

Amylin antibody

Catalog Number: 108015

Product name

Amylin antibody

Specificity

Human; other species not tested.

Antibody description

Amylin Rabbit Polyclonal antibody. Positive IHC detected in human pancreas tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000415).
Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:108015(Amylin Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:108015(Amylin Antibody) at dilution of 1:50 (under 40x lens)