

Anti-FADD antibody



Catalog Number: 100415

Product name

Anti-FADD antibody

Specificity

Human FADD

Antibody description

Rabbit polyclonal to FADD

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human FADD, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 10-20 µg/ml


IP: 0.5-2 uL/mg of lysate

Validations



Items	Lanes	A
Sample (whole cell lysate)		Jurkat
Sample Volume (µg/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		10-20 µg/ml
Secondary Antibody		Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.
		Developed using Odyssey imaging system.
Explanation		Predicted band size : 23 kDa Observed band size : 27 kDa Additional bands at : 37 kDa [We are unsure as to the identity of these extra bands.]

FADD Antibody, Rabbit PAb, Antigen Affinity



Items	Lanes	A	B
Sample (whole cell lysate)		HeLa	Jurkat
Sample quantity		0.5 mg	
IP antibody quantity		2 µg	
Immunomagnetic beads Protein G		60 µg	
Gel		13% SDS-PAGE reducing gel	
Primary antibody		KLJ1-SMCC-sinoA909 antibody at 10 µg/ml	
Secondary antibody		Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:5000 dilution	

Purified, Western blot

FADD Antibody, Rabbit PAb, Antigen Affinity

Purified, Immunoprecipitation