

LYVE1 antibody



Catalog Number: 112415

Product name

LYVE1 antibody

Specificity

Human; other species not tested.

Antibody description

LYVE1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of LYVE1 recombinant protein (Accession Number: NM_006691). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112415(LYVE1 antibody) at dilution of 1:500