

Anti-CCL4L1 antibody

Catalog Number: 176415

Product name

Anti-CCL4L1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CCL4L1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human ccl4. aa range:41-90

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of HL-60, 22RV-1 cells using MIP-1 β Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000