

Anti-CCL4L1 antibody



Catalog Number: 176415

Product name

Polyclonal

Anti-CCL4L1 antibody

Ig Type

Specificity

Rabbit IgG

Human

Applications

Antibody description

WB

Rabbit polyclonal antibody to CCL4L1

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human ccl4. aa range:41-90

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1.; Western Blot analysis of HL-60, 22RV-1 cells using MIP-1β Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000