

# TRPS1 antibody

Catalog Number: 116415

## Product name

TRPS1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TRPS1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, HEK-293 cells, HeLa cells. Positive IP detected in mouse lung tissue. Observed molecular weight by Western-blot: 150 kDa

## Preparation

This antibody was obtained by immunization of TRPS1 recombinant protein (Accession Number: BC125020). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:116415 (TRPS1 antibody) at dilution of 1:1000



IP Result of anti-TRPS1 (IP: Catalog No: 116415, 4ug; Detection: Catalog No: 116415 1:1000) with mouse lung tissue lysate 3600ug.