

Anti-Cyclin A1 Rabbit antibody



Catalog Number: 167415

Product name

Anti-Cyclin A1 Rabbit antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Cyclin A1 Rabbit polyclonal antibody

Preparation

Antigen: The antiserum was produced against synthesized peptide derived from human Cyclin A. AA range:221-270

Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, ELISA

Dilutions


WB: 1/500 - 1/2000

IHC: 1/100 - 1/300

IF: 1/200 - 1/1000

ELISA: 1/10000

Validations



Western Blot analysis of various cells using Cyclin A Polyclonal Antibody