

MAPK13 antibody



Catalog Number: 107415

Product name

MAPK13 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MAPK13 Mouse Monoclonal antibody. Positive WB detected in HepG2 cells, HeLa cells, MCF-7 cells. Observed molecular weight by Western-blot: 42kd

Preparation

This antibody was obtained by immunization of MAPK13 recombinant protein (Accession Number: NM_002754). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107415(MAPK13 antibody) at dilution of 1:1000