

Human ATG4A/APG4A/AUTL2 (His tag) recombinant protein



Catalog Number: 518415

General Information

Gene Name Synonym

Cysteine Protease ATG4A, AUT-Like 2 Cysteine Endopeptidase, Autophagin-2, Autophagy-Related Cysteine Endopeptidase 2, Autophagy-Related Protein 4 Homolog A, hAPG4A, ATG4A, APG4A, AUTL2

Protein Construction

A DNA sequence encoding the Human ATG4A/APG4A/AUTL2 Met1-Val398 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ATG4A/APG4A/AUTL2 (His tag)
recombinant protein