

GNAZ Polyclonal Antibody



Catalog Number: 160115

Product name

GNAZ Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to GNAZ

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-355 of human GNAZ (NP_002064.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - GNAZ Polyclonal Antibody

Western blot analysis of extracts of mouse brain, using GNAZ antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.