

# GPR153 antibody

Catalog Number: 111115

## Product name

GPR153 antibody

## Specificity

Human; other species not tested.

## Antibody description

GPR153 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive WB detected in Jurkat cells, K-562 cells. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of GPR153 recombinant protein (Accession Number: XM\_011541434). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:50-1:500


## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111115(GPR153 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:111115(GPR153 Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:111115(GPR153 Antibody) at dilution of 1:100 (under 40x lens).