

Anti-NEUROD1 antibody



Catalog Number: 106115

Product name

Anti-NEUROD1 antibody

Specificity

Human NEUROD1

Antibody description

Rabbit polyclonal to NEUROD1

Preparation

Produced in rabbits immunized with A synthetic peptide corresponding to the N-terminus of the human NeuroD1, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 0.1-2 µg/mL

Validations



NeuroD1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human NEUROD1 in human hepatoma with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded



sections).

NeuroD1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human NEUROD1 in human endometrium with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded sections).