

AP1B1 antibody

Catalog Number: 108115

Product name

AP1B1 antibody

IP: 1:200-1:2000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AP1B1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HEK-293 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 101 kDa

Preparation

This antibody was obtained by immunization of AP1B1 recombinant protein (Accession Number: NM_001166019). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108115(AP1B1 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using AP1B1 antibody Catalog No:108115 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-AP1B1 (IP:Catalog No:108115,

AP1B1 antibody



Catalog Number: 108115

4ug; Detection:Catalog No:108115 1:500) with HEK-293 cells lysate 2000ug.