

Mouse Ccl2/Je/Mcp1/Scya2 (His tag) recombinant protein



Catalog Number: 518115

General Information

Gene Name Synonym

C-C motif chemokine 2; Monocyte chemoattractant protein 1; Monocyte chemotactic protein 1; MCP-1; Platelet-derived growth factor-inducible protein JE; Small-inducible cytokine A2; Ccl2; Je; Mcp1; Scya2

Protein Construction

A DNA sequence encoding the Mouse Ccl2/Je/Mcp1/Scya2 Gln24-Asn148 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 500mM NaCl, 10% glycerol, pH7.4 .

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Ccl2/Je/Mcp1/Scya2 (His tag) recombinant protein