

HNRNPUL2-Specific antibody



Catalog Number: 111515

Product name

HNRNPUL2-Specific antibody

Specificity

Human; other species not tested.

Antibody description

HNRNPUL2-Specific Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in MCF-7 cells. Observed molecular weight by Western-blot: 100-110 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001079559). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:50-1:500

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:111515(HNRNPUL2-Specific Antibody) at dilution of 1:300



Immunofluorescent analysis of (10% Formaldehyde) fixed MCF-7 cells using Catalog No:111515(HNRNPUL2-Specific Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)