

# OTUD5 antibody



Catalog Number: 113515

## Product name

OTUD5 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

OTUD5 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, HeLa cells. Positive IF detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 61 kDa

## Preparation

This antibody was obtained by immunization of OTUD5 recombinant protein (Accession Number: BC028225). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

## Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113515(OTUD5 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (10% Formaldehyde) fixed NIH/3T3 cells using Catalog No:113515(OTUD5 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)