

# Anti-Hsc70 Mouse antibody



Catalog Number: 168515

## Product name

Anti-Hsc70 Mouse antibody

## Specificity

Human

## Antibody description

Hsc70 Mouse monoclonal antibody

## Preparation

Antigen: Recombinant Protein of Hsc 70.

## Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide 50% glycerol.

## Storage

Store at -20°C

## Clonality

Monoclonal

## Ig Type

Mouse IgG


## Applications

WB

## Dilutions

1:1000

## Validations



Western blot detection of Hsc70 in HeLa,293 cell lysates using Hsc70 Mouse mAb(1:1000 diluted).Predicted band size:70KDa.Observed band size:70KDa.