

AVPR1B antibody

Catalog Number: 108515

Product name

AVPR1B antibody

Specificity

Human; other species not tested.

Antibody description

AVPR1B Rabbit Polyclonal antibody. Positive IHC detected in human pancreas tissue, human kidney tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000707).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:117076(AVPR1B Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:117076(AVPR1B Antibody) at dilution of 1:100 (under 40x lens).