

KATNAL1 antibody



Catalog Number: 111915

Product name

KATNAL1 antibody

Immunogen

[Human KATNAL1 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

KATNAL1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 50-55 kDa

Preparation

This antibody was obtained by immunization of KATNAL1 recombinant protein (Accession Number: XM_024449423). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111915(KATNAL1 Antibody) at dilution of



1:300

IP Result of anti-KATNAL1 (IP:Catalog No:111915, 3ug; Detection:Catalog No:111915 1:500) with HeLa cells lysate 3000ug.