

Anti-Phospho-Rsk 2/MAPKAP Kinase 1b (Ser227) Rabbit antibody



Catalog Number: 166915

Product name

Polyclonal

Anti-Phospho-Rsk 2/MAPKAP Kinase 1b (Ser227)
Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IP

Antibody description

Phospho-Rsk 2/MAPKAP Kinase 1b (Ser227)
Rabbit polyclonal antibody

Dilutions


WB: 1/1000

IP: 1/20

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser227 of human Rsk 2/MAPKAP Kinase 1b

Validations



Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Western blot detection of Phospho-Rsk 2/MAPKAP Kinase 1b (Ser227) in 3T3,HeLa cell lysates using Phospho-Rsk 2/MAPKAP Kinase 1b (Ser227) Rabbit pAb(1:1000 diluted).Predicted band size:84kDa.Observed band size:84kDa.