

# Anti-HDAC9 Rabbit antibody



Catalog Number: 167915

## Product name

Anti-HDAC9 Rabbit antibody

## Specificity

Human, Rat

## Antibody description

HDAC9 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human HDAC9

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/2000-1/10000

IHC: 1/50

ICC/IF: 1/50

## Validations



Western blot detection of HDAC9 in A549,HL-60 cell lysates using HDAC9 Rabbit mAb(1:1000 diluted).Predicted band size:111kDa.Observed band size:150kDa.