

Anti-Bax Rabbit antibody

Catalog Number: 168915

Product name

Rabbit IgG

Anti-Bax Rabbit antibody

Applications

Specificity

WB, IHC-P, IP

Human, Mouse, Rat

Dilutions

Antibody description

1:1000

Bax Rabbit monoclonal antibody

1:20-1:500

Preparation

1:20

Antigen: A synthetic peptide of human Bax

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage


Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Validations



Western blot detection of Bax in C6,3T3,HeLa cell lysates using Bax Rabbit mAb(1:1000 diluted).Predicted band size:21kDa.Observed band size:20kDa.