

# Human CCNH (His tag) recombinant protein



Catalog Number: 518915

## General Information

### Gene Name Synonym

Cyclin-H,CCNH,MO15-associated protein,p34,p37

### Protein Construction

A DNA sequence encoding the Human CCNH Met1-Leu323 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

## Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

## Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM Tris, 100mM NaCl, 2mM EDTA, 2mM DTT, 30% Glycerol, pH8.5.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

## SDS-PAGE



Human CCNH (His tag) recombinant protein