

TPRKB antibody



Catalog Number: 116215

Product name

TPRKB antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

TPRKB Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Observed molecular weight by Western-blot: 20kd

Preparation

This antibody was obtained by immunization of TPRKB recombinant protein (Accession Number: NM_001330389). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116215(TPRKB antibody) at dilution of 1:300