

# Anti-FEN1 (1E7) Mouse antibody



Catalog Number: 168215

## Product name

Anti-FEN1 (1E7) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

FEN1 (1E7) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human FEN-1 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB, IP


## Dilutions

WB: 1/1000


## Validations



Western blot detection of FEN-1 in HeLa, HeLa NE, Jurkat, K562 and 3T3 cell lysates using FEN-1 mouse mAb (1:1000 diluted). Predicted band size: 45 KDa. Observed band size: 45 KDa.



Western blot detection of FEN-1 in HeLa, Jurkat, 3T3, C6, CEM and Ramos cell lysates using FEN-1 mouse mAb (1:1000 diluted). Predicted band size: 45 KDa. Observed



band size: 45 KDa.

Immunoprecipitation analysis of HeLa cell lysates using FEN-1 mouse mAb.