

# CAV2 Monoclonal Antibody



Catalog Number: 160615

## Product name

CAV2 Monoclonal Antibody

## Specificity

Human

## Antibody description

Monoclonal antibody to CAV2

## Preparation

Antigen: Recombinant protein of human CAV2

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF, IP

## Dilutions


WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200


IP 1:20 - 1:50

## Validations




Western blot - CAV2 Monoclonal Antibody

Western blot analysis of extracts of HUVEC cells, using CAV2 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - CAV2 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded human cervix using CAV2 antibody at dilution of 1:100 (40x lens).



# CAV2 Monoclonal Antibody




Catalog Number: 160615

---

Immunohistochemistry - CAV2 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded human uterus using CAV2 antibody at dilution of 1:100 (40x lens).



Immunofluorescence - CAV2 Monoclonal Antibody

Immunofluorescence analysis of HeLa cells using CAV2 antibody .