

Cynomolgus Interferon alpha-B / IFNA8 (His Tag) recombinant protein



Catalog Number: 500615

General Information

Protein Construction

A DNA sequence encoding the cynomolgus IFNA8 (G7PS89) (Met1-Glu189) was expressed with a polyhistidine tag at the C-terminus.

Organism

Cynomolgus

Expression Host

Human Cells

QC Testing

Activity

Measured in antiviral assays using WISH human amnion cells infected with vesicular stomatitisvirus (VSV).

The ED₅₀ for this effect is 0.2-1.6 pg/mL.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 33

Molecular Mass

The recombinant cynomolgus IFNA8 comprises 177 amino acids and has a calculated molecular mass of 21 KDa. The apparent molecular mass of it is approximately 33 KDa respectively in SDS-PAGE.

Formulation

Lyophilized from sterile PBS, PH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Cynomolgus Interferon alpha-B / IFNA8 Protein (His Tag) SDS-PAGE