

IFI30 antibody

Catalog Number: 111615

Product name

IFI30 antibody

Specificity

Human; other species not tested.

Antibody description

IFI30 Rabbit Polyclonal antibody. Positive WB detected in human colon tissue. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of IFI30 recombinant protein (Accession Number: NM_006332). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:111615(IFI30 antibody) at dilution of 1:400