

MIDN antibody



Catalog Number: 112615

Product name

MIDN antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MIDN Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of MIDN recombinant protein (Accession Number: NM_177401). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:112615(MIDN Antibody) at dilution of 1:300