

Anti-Neuropilin-2/NRP2 antibody



Catalog Number: 104615

Product name

Anti-Neuropilin-2/NRP2 antibody

Immunogen

[Human Neuropilin-2/NRP2 \(His Tag\) recombinant protein](#)

Specificity

Human Neuropilin 2 / NRP2

Antibody description

Mouse monoclonal to Neuropilin-2/NRP2

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human Neuropilin 2 / NRP2 (rh Neuropilin 2 / NRP2; NP_003863.2; Met 1-Tyr 855) and conjugated with APC under optimum conditions, the unreacted APC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.1% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze !

Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



Neuropilin 2 / NRP2 Antibody (APC), Mouse MAb, Flow cytometric analysis

Profile of anti-NRP2 reactivity on HUVEC cells analyzed by flow cytometry.