

# CACNA1B Polyclonal Antibody



Catalog Number: 163315

## Product name

CACNA1B Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to CACNA1B

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 2200 to the C-terminus of human CACNA1B (NP\_000709.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - CACNA1B Polyclonal Antibody

Western blot analysis of extracts of U-87MG cells, using CACNA1B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3min.