

# NOX5 antibody



Catalog Number: 113315

## Product name

NOX5 antibody

## Immunogen

[Human NOX5 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

NOX5 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, Raji cells. Observed molecular weight by Western-blot: 60 kDa

## Preparation

This antibody was obtained by immunization of NOX5 recombinant protein (Accession Number: BC125098). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113315(NOX5 Antibody) at dilution of 1:1000