

TRIM39 antibody



Catalog Number: 116315

Product name

TRIM39 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

TRIM39 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 45kd,56-60kd

Preparation

This antibody was obtained by immunization of TRIM39 recombinant protein (Accession Number: NM_001369523). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116315(TRIM39 antibody) at dilution of 1:300