

Anti-CTNNA1 Rabbit antibody



Catalog Number: 167315

Product name

Anti-CTNNA1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

CTNNA1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human alpha 1 Catenin

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/5000


IHC: 1/20

ICC/IF: 1/50

FC: 1/20

IP: 1/20

Validations



Western blot detection of alpha 1 Catenin in Rat Brain, C6, 3T3, HeLa cell lysates using alpha 1 Catenin Rabbit pAb(1:1000 diluted). Predicted band size: 100kDa. Observed band size: 100kDa.