

# ATP6V1G2 antibody

Catalog Number: 108315

## Product name

ATP6V1G2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ATP6V1G2 Rabbit Polyclonal antibody. Positive IHC detected in mouse brain tissue. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 14 kDa

## Preparation

This antibody was obtained by immunization of ATP6V1G2 recombinant protein (Accession Number: NM\_130463). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108315(ATP6V1G2 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:108315(ATP6V1G2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:108315(ATP6V1G2 Antibody) at dilution of 1:50 (under 40x lens)