

TYR Polyclonal Antibody



Catalog Number: 161715

Product name

TYR Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to TYR

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-340 of human TYR (NP_000363.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF


Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200


Validations

Rat skin



Western blot - TYR Polyclonal Antibody

Western blot analysis of extracts of rat skin, using TYR antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time:



Immunofluorescence - TYR Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using TYR antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.