

# Cystatin C antibody



Catalog Number: 109715

## Product name

Cystatin C antibody

IHC: 1:20-1:200

IF: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Validations

## Antibody description

Cystatin C Rabbit Polyclonal antibody. Positive IHC detected in human ovary tissue, human brain tissue. Positive IF detected in HeLa cells. Positive WB detected in human brain tissue, fetal human brain tissue. Observed molecular weight by Western-blot: 16 kDa




## Preparation

This antibody was obtained by immunization of Cystatin C recombinant protein (Accession Number: NM\_000099). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109715(CST3 antibody) at dilution of 1:400

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human ovary tissue slide using Catalog No:109715(Cystatin C Antibody) at dilution of 1:200 (under 10x lens).

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IF, IHC

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunohistochemistry of paraffin-embedded


# Cystatin C antibody



Catalog Number: 109715

---

human ovary tissue slide using Catalog No:109715(Cystatin C Antibody) at dilution of 1:200 (under 40x lens).



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:109715(Cystatin C Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)