

# GLIS1 antibody



Catalog Number: 111025

## Product name

GLIS1 antibody

## Specificity

Human; other species not tested.

## Antibody description

GLIS1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 66 kDa

## Preparation

This antibody was obtained by immunization of GLIS1 recombinant protein (Accession Number: NM\_147193). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111025(GLIS1 Antibody) at dilution of 1:300