

# Anti-HOXA5 antibody



Catalog Number: 176025

## Product name

Anti-HOXA5 antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to HOXA5

## Preparation

This antigen of this antibody was synthesized peptide derived from the n-terminal region of human hoxa5.

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1;; Western Blot analysis of 293T cells using HoxA5 Polyclonal Antibody diluted at 1: 500 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).