

TGFBRAP1 antibody



Catalog Number: 116025

Product name

TGFBRAP1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human; other species not tested.

Validations

Antibody description

TGFBRAP1 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Positive IF detected in HepG2 cells. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 80 kDa




Preparation

This antibody was obtained by immunization of TGFBRAP1 recombinant protein (Accession Number: XM_047446423). Purification method: Antigen affinity purified.

COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:116025(TGFBRAP1 Antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:116025(TGFBRAP1 Antibody) at dilution of 1:50 (under 40x lens)



TGFBRAP1 antibody



Catalog Number: 116025

Immunofluorescent analysis of HepG2 cells using
Catalog No:116025(TGFBRAP1 Antibody) at

dilution of 1:25 and Alexa Fluor 488-conjugated
AffiniPure Goat Anti-Rabbit IgG(H+L)