

Anti-NADPH Oxidase 4 Rabbit antibody



Catalog Number: 167025

Product name

Anti-NADPH Oxidase 4 Rabbit antibody

ICC/IF: 1/50-1/200

IP: 1/50

Specificity

Human, Mouse, Rat


Validations

Antibody description

NADPH Oxidase 4 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human NOX4




Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Western blot analysis of NOX4 expression in JAR

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.



Clonality

Polyclonal

cell lysate.


Ig Type

Rabbit IgG

Immunohistochemical analysis of paraffin-embedded rat kidney, using NOX4 Antibody .

Applications

WB, IHC-P, ICC/IF, IP



Dilutions

WB: 1/1000-1/2000

IHC: 1/50-1/200

Immunofluorescent analysis of HeLa cells, using NOX4 Antibody .