

# BSAP,PAX5 antibody



Catalog Number: 107025

## Product name

BSAP,PAX5 antibody

## Specificity

Human; other species not tested.

## Antibody description

BSAP,PAX5 Mouse Monoclonal antibody. Positive WB detected in Ramos cells, Daudi cells, Raji cells. Observed molecular weight by Western-blot: 42 kDa

## Preparation

This antibody was obtained by immunization of BSAP,PAX5 recombinant protein (Accession Number: NM\_016734). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Ramos cells were subjected to SDS PAGE followed by western blot with Catalog No:107025(BSAP,PAX5 Antibody) at dilution of 1:1500