

Human LYPLA1/APT1/LPL1 (His tag) recombinant protein



Catalog Number: 518025

General Information

Gene Name Synonym

Acyl-Protein Thioesterase 1, APT-1, hAPT1, Lysophospholipase 1, Lysophospholipase I, LPL-I, LysoPLA I, LYPLA1, APT1, LPL1

Protein Construction

A DNA sequence encoding the Human LYPLA1/APT1/LPL1 Met1-Asp230 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human LYPLA1/APT1/LPL1 (His tag) recombinant protein