

Human FIBP recombinant protein



Catalog Number: 519025

General Information

Gene Name Synonym

Acidic Fibroblast Growth Factor Intracellular-Binding Protein, aFGF Intracellular-Binding Protein, FGF-1 Intracellular-Binding Protein, FIBP

Protein Construction

A DNA sequence encoding the Human FIBP Met1-Asp357 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mM TrisHCl pH8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human FIBP recombinant protein