

ORM1 antibody



Catalog Number: 113425

Product name

ORM1 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200


Specificity

Human, Mouse; other species not tested.

Validations

Antibody description

ORM1 Rabbit Polyclonal antibody. Positive IP detected in human plasma tissue. Positive WB detected in mouse liver tissue, mouse bladder tissue. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 41-47kd




Preparation

This antibody was obtained by immunization of ORM1 recombinant protein (Accession Number: NM_000607). Purification method: Antigen affinity purified.

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113425(ORM1 antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human liver slide using Catalog No:113425(ORM1 Antibody) at dilution of 1:50

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

Immunohistochemistry of paraffin-embedded human liver slide using Catalog No:113425(ORM1

ORM1 antibody

Catalog Number: 113425

Antibody) at dilution of 1:50



IP Result of anti-ORM1 (IP:Catalog No:113425, 3ug; Detection:Catalog No:113425 1:300) with human plasma tissue lysate 100ug.