

Anti-Cyclin E1 Rabbit antibody



Catalog Number: 167425

Product name

Anti-Cyclin E1 Rabbit antibody

Specificity

Human

Antibody description

Cyclin E1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Cyclin E1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/50

Validations



Western blot detection of Cyclin E1 in K562, HeLa cell lysates using Cyclin E1 Rabbit pAb(1:1000 diluted). Predicted band size: 47kDa. Observed band size: 47kDa.