

# CNKSR1 antibody

Catalog Number: 109425

## Product name

CNKSR1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

CNKSR1 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 80kd

## Preparation

This antibody was obtained by immunization of CNKSR1 recombinant protein (Accession Number: NM\_001297647). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109425(CNKSR1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:109425(CNKSR1 antibody) at dilution of 1:100 (under 40x lens)