

# Anti-RAB10 Rabbit antibody



Catalog Number: 166825

## Product name

Anti-RAB10 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

RAB10 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human RAB10

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/1000

ICC/IF: 1/100

FC: 1/100

IP: 1/20

## Validations



Western blot detection of RAB10 in C6,3T3,Hela cell lysates using RAB10 Rabbit pAb(1:1000 diluted).Predicted band size:23kDa.Observed band size:23kDa.