

Anti-Mras Rabbit antibody



Catalog Number: 167825

Product name

Anti-Mras Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Mras Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Mras

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/1000

IP: 1/20

Validations



Western blot detection of Mras in HeLa,Rat Brain lysates using Mras Rabbit mAb(1:1000 diluted).Predicted band size:24kDa.Observed band size:24kDa.