

# Anti-PCBP2/hnRNP E2 Rabbit antibody



Catalog Number: 168825

## Product name

Rabbit IgG

Anti-PCBP2/hnRNP E2 Rabbit antibody

## Applications

## Specificity

WB, ICC/IF, FC, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

1:1000

PCBP2/hnRNP E2 Rabbit monoclonal antibody

1:100

## Preparation

1:50


Antigen: A synthetic peptide of human PCBP2

1:20

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Validations



## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of PCBP2/hnRNP E2 in C6,3T3 cell lysates using PCBP2/hnRNP E2 Rabbit mAb(1:1000 diluted).Predicted band size:39kDa.Observed band size: 35-45kDa.

## Clonality

Monoclonal

## Ig Type